

# SANI-CLOTH®

## Disinfection of Computer Keyboards



**What's growing  
on your keyboard?**

### Recommendations for disinfection of computer keyboards<sup>1,2</sup>:

In the office, patient testing area, nursing unit or at the patient bedside, the computer keyboard can serve as a reservoir for pathogenic microorganisms.

- Transmission of pathogens present in environments of patients infected or colonized with MRSA, VRE and *Clostridia difficile* are of special concern
- Computer keyboards can serve as reservoirs for cross transmission of microorganisms such as bacteria, viruses and fungi
- Contaminated computer keyboards can serve as a mechanism for contaminating hands of healthcare workers
- Studies have shown microbial contamination of computer keyboards is prevalent<sup>1</sup>
- Computer keyboards can be successfully decontaminated using Sani-Cloth®
- Sani-Cloth® demonstrated excellent sustained activity against VRE and antimicrobial activity was maintained over the entire 48-hour test period<sup>1</sup>

- Always check with equipment manufacturer for recommended disinfectant chemicals
- Do not touch keyboards with contaminated hands
- Disinfect computer keyboards daily and when visibly soiled
- Computers should remain in rooms of patients in isolation precautions
- Mobile computers should be disinfected prior to entry and upon exit of patient room/area















<sup>1</sup>Rutala et al. "Bacterial Contamination of Keyboards: Efficacy and Functional Impact of Disinfectants." ICHE Vol. 27, No. 4. April, 2006. pp. 372-377.





<sup>2</sup>Neely et. al. "Computer equipment used in patient care within a multihospital system and recommendations for cleaning and disinfection." APIC Practice Forum. AJIC. Vol. 93, No. 4., May, 2005. pp. 233-237.



# SANI-CLOTH®

## Disinfection of Computer Keyboards

Description	Reorder No.	Wipe Size	Case Pack	Case Weight	Case Dimension L x W x H	Case Cube	Pallet TI/HI
<b>WIPES</b>							
 Sani-Cloth® Plus Large	Q89072	6 x 6.75 in 15.25 x 17.14 cm	160/can 12 can/cs	26 lbs 11.793 kg	15.8125 x 10.8125 x 14.50 in 40.16 x 27.46 x 36.80 cm	1.43 cu ft 0.046 cu m	10/3
 Sani-Cloth® Plus Extra Large	Q85084	7.5 x 15 in 19.05 x 38.10 cm	65/can 6 can/cs	13.80 lbs 6.259 kg	15.6875 x 10.875 x 10.25 in 39.84 x 27.62 x 26.03 cm	1.012 cu ft 0.0287 cu m	10/4
 Sani-Cloth® Plus Individual Large	H54082	5 x 8 in 12.70 x 20.32 cm	50/box 10 box/cs	7.72 lbs 3.501 kg	13.125 x 7 x 9.50 in 26.03 x 17.78 x 24.13 cm	0.5051 cu ft 0.0143 cu m	19/4
 Sani-Cloth® Plus Individual Extra Large	U36495	11.5 x 11.75 in 29.21 x 29.84 cm	50/box 3 box/cs	7.16 lbs 3.247 kg	13.9375 x 4.8125 x 12.875 in 35.40 x 12.22 x 32.70 cm	0.4998 cu ft 0.0142 cu m	24/3
 Super Sani-Cloth® Large	Q55172	6 x 6.75 in 15.25 x 17.14 cm	160/can 12 can/cs	22.60 lbs 10.273 kg	15.8125 x 10.8125 x 14.50 in 40.16 x 27.46 x 36.80 cm	1.43 cu ft 0.046 cu m	10/3
 Super Sani-Cloth® Extra Large	Q86984	7.5 x 15 in 19.05 x 38.10 cm	65/can 6 can/cs	14.30 lbs 6.486 kg	15.6875 x 10.875 x 10.25 in 39.84 x 27.62 x 26.03 cm	1.012 cu ft 0.0287 cu m	10/4
 Super Sani-Cloth® Individual Large	H04082	5 x 8 in 12.70 x 20.32 cm	50/box 10 box/cs	7.72 lbs 3.501 kg	13.125 x 7 x 9.50 in 26.03 x 17.78 x 24.13 cm	0.5051 cu ft 0.0143 cu m	19/5
 Super Sani-Cloth® Individual Extra Large	U87295	11.50 x 11.75 in 29.21 x 29.84 cm	50/box 3 box/cs	8.00 lbs 3.628 kg	13.9375 x 4.8125 x 12.875 in 35.40 x 12.22 x 32.70 cm	0.4998 cu ft 0.0142 cu m	24/3
 Sani-Cloth® HB Large	Q08472	6 x 6.75 in 15.25 x 17.14 cm	160/can 12 can/cs	24.75 lbs 11.226 kg	15.8125 x 10.8125 x 14.50 in 40.16 x 27.46 x 36.8 cm	1.43 cu ft 0.046 cu m	10/3
 Sani-Cloth® HB Extra Large	Q85484	7.5 x 15 in 19.05 x 38.10 cm	65/can 6 can/cs	14.30 lbs 6.486 kg	15.6875 x 10.875 x 10.25 in 39.84 x 27.62 x 26.03 cm	1.012 cu ft 0.0287 cu m	10/4
 Sani-Cloth® HB Individual Large	H21082	5 x 8 in 12.70 x 20.32 cm	50/box 10 box/cs	7.72 lbs 3.501 kg	13.125 x 7 x 9.50 in 26.03 x 17.78 x 24.13 cm	0.5051 cu ft 0.0143 cu m	19/5
 Sani-Cloth® HB Individual Extra Large	U22495	11.50 x 11.75 in 29.21 x 29.84 cm	50/box 3 box/cs	7.90 lbs 3.583 kg	13.9375 x 4.8125 x 12.875 in 35.40 x 12.22 x 32.70 cm	0.4998 cu ft 0.0142 cu m	24/3

Description	Reorder No.	Use with Can Size	Case Weight	Case Dimension L x W x H	Case Cube	Pallet TI/HI
<b>SANI-BRACKET</b>						
 Sani-Bracket® 3-in-1 (Sani-Cloth®) Labeled: Not for Use on Skin	P010900	Large	19.50 lbs	29.75 x 12.75 x 9.25 in 75.5 x 32.5 x 23.5 cm	2.03 cu ft	4/4
 Sani-Bracket® 3-in-1 (Sani-Cloth®) Labeled: Not for Use on Skin	P010901	Extra Large	1.75 lbs	6 x 6 x 8.25 in 15.25 x 15.25 x 21 cm	0.17 cu ft	42/4
 Sani-Bracket® 3-in-1 Pole Accessory Kit	P010600	P010900 P010901	1.65 lbs	5.5 x 5.25 x 3 in 13.97 x 13.335 x 7.62 cm	0.065 cu ft	63/14
 Sani-Bracket® for X-Large Box (Sani-Cloth®)	P44210	U87295 U22495	0.60 lbs	5.25 x 5.25 x 5.25 in 13.5 x 13.5 x 13.5 cm	0.07 cu ft	72/8

**GSA** U.S. General Services Administration  
Sani-System® Federal GSA# GS-07F-0412J.

Professional Disposables International, Inc.  
Two Nice-Pak Park  
Orangeburg, NY 10962-1376  
ISO 9001:2008 Certified  
Fax: (845) 398-5347  
1-800-999-6423 in USA  
[www.pdipi.com](http://www.pdipi.com)

©2010 Professional Disposables International, Inc.  
Sani-System®, Sani-Bracket® and  
Sani-Cloth® are registered trademarks of  
Professional Disposables International, Inc.

