

# SANI-CLOTH® BRAND

## Germicidal Disposable Wipes



### TECHNICAL SUMMARY AND DATA

Features:	SUPER	AF	PLUS
Recommended contact time	2 minutes	3 minutes	5 minutes
Alcohol content	55.0 %	Alcohol Free*	14.85%
Meets CDC Disinfection and Sterilization Guidelines <sup>1</sup>	●	●	●
EPA registered hospital grade disinfectant	●	●	●
Quick pull canisters - Large wipe (160 count)	6.0" x 6.75"	6.0" x 6.75"	6.0" x 6.75"
Quick pull canisters - Extra large wipe (65 count)	7.5" x 15.0"	7.5" x 15.0"	7.5" x 15.0"
Individual packets - Large wipe (50 count)	5.0" x 8.0"	5.0" x 8.0"	5.0" x 8.0"
Individual packets - Extra large wipe (50 count)	11.5" x 11.75"	11.5" x 11.75"	11.5" x 11.75"
Compliance accessories to meet The Joint Commission <sup>2</sup> & CMS requirements <sup>3</sup>	●	●	●
Tested Effective Against:			
<i>Acinetobacter baumannii</i>	●	●	
Adenovirus	●		
Avian Influenza A Virus(H5N1)		●	
<i>Burkholderia cepacia</i>	●		
<i>Campylobacter jejuni</i>	●	●	●
<i>Candida albicans</i>	●	●	
Community Acquired Methicillin Resistant <i>Staphylococcus aureus</i> (CA-MRSA) Genotype USA300		●	
Community Acquired Methicillin Resistant <i>Staphylococcus aureus</i> (CA-MRSA) Genotype USA400		●	
ESBL Producing <i>Escherichia coli</i> (E. coli)	●	●	
<i>Escherichia coli</i> (E. coli)	●	●	●
<i>Escherichia coli</i> (E. coli) O157:H7	●	●	●
Hepatitis B Virus	●	●	●
Hepatitis C Virus	●	●	●
Herpes Simplex Virus, Type 2	●	●	●
HIV-1 (Aids Virus)	●	●	●
Human Coronavirus	●	●	
Pandemic 2009 H1N1 Influenza A Virus	●	●	●
Influenza A2 / Hong Kong Virus	●	●	●
<i>Klebsiella pneumoniae</i>	●	●	
<i>Listeria monocytogenes</i>		●	
Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA)	●	●	●
<i>Pseudomonas aeruginosa</i>	●	●	●
Respiratory Syncytial Virus	●	●	●
Rhinovirus	●		
Rotavirus	●		
<i>Salmonella enterica</i> (formerly <i>Salmonella choleraesuis</i> )	●	●	●
<i>Staphylococcus aureus</i>	●	●	●
<i>Streptococcus pyogenes</i>		●	
TB ( <i>Mycobacterium bovis</i> )	●		●
Vaccinia Virus	●		
Vancomycin Resistant <i>Enterococcus faecalis</i> (VRE)	●	●	●
<b>Total Microorganisms Tested</b>	<b>26</b>	<b>25</b>	<b>16</b>

\* contains less than 0.016% alcohol <sup>1</sup> Guideline for Disinfection and Sterilization in Healthcare Facilities, Centers for Disease Control and Prevention, 2008 <sup>2</sup> Infection Prevention and Control (IPC) Standard IC.02.02.01 and Elements of Performance (EP) 1 and 2, The Joint Commission, October 2009 <sup>3</sup> Memorandum: Point of Care Devices and Infection Control in Nursing Homes, Center for Medicaid, CHIP, and Survey & Certification/Survey & Certification Group, August 27, 2010

© 2011 Professional Disposables International, Inc.  
Sani-Cloth® is a registered trademark of Professional Disposables International, Inc.

Professional Disposables International, Inc.

ISO 9001:2008 Certified

1-800-999-6423 in USA

www.pdipdi.com

PDI12104095

