

SANI-CLOTH® BRAND

Germicidal Disposable Wipes

TECHNICAL SUMMARY & DATA



Features:	PLUS	SUPER	HB
Intermediate-level disinfectant	●	●	
Quick pull canisters Large wipe 6 x 6.75 in. (160 count)	●	●	●
Quick pull canisters X-Large wipe 8 x 14 in. (65 count)	●	●	●
Pocket-sized packets	●	●	●
PPM (Parts Per Million) quaternary ammonium chlorides	2,500	5,000	1,400
Alcohol content	LOW: 14.85%	HIGH: 55%	LESS THAN 0.15%
Recommended contact time (*varies by microorganism – see Technical Data Bulletin)	5 min. or less*	2 min. or less*	10 min. or less*
Sani-Bracket® 3-in-1 Wall/Mobile Equipment Attachment System	●	●	●

Tested Effective Against:

TB (Mycobacterium bovis)	●	●	
Hepatitis B Virus	●	●	●
Hepatitis C Virus	●	●	
HIV-1	●	●	●
Herpes Simplex Type 2	●	●	●
Respiratory Syncytial Virus	●	●	
Methicillin Resistant Staphylococcus aureus (MRSA)	●	●	●
Vancomycin Resistant Enterococcus (VRE)	●	●	●
Avian Influenza A (H1N1)		●	
Influenza A2	●	●	●
Escherichia coli (E. coli)	●	●	●
Escherichia coli O157:H7 (E. coli)	●	●	●
ESBL Producing Escherichia coli (E. coli)		●	
Pseudomonas aeruginosa	●	●	●
Salmonella enterica (formerly Salmonella choleraesuis)	●	●	●
Staphylococcus aureus	●	●	●
Klebsiella pneumoniae		●	●
Vaccinia		●	●
Candida albicans		●	●
Rhinovirus		●	
Human Coronavirus		●	
Rotavirus		●	
Adenovirus		●	●
Burkholderia cepacia		●	●
Acinetobacter baumannii		●	
Campylobacter jejuni		●	●
Total Microorganisms Tested	14	26	100+
*Refer to Technical Data Bulletin for complete listing	●	●	●

Professional Disposables International, Inc.
 ISO 9001:2000 Certified
 1-800-999-6423 in USA
 www.pdiipi.com

